THC

EasyTherm

DE	<i>Bedienungsanleitung – EasyTherm Ohrthermometer</i>	A - 1
EN	<i>Operating manual – EasyTherm Ear Thermometer</i>	B - 1
FR	<i>Manuel d'utilisation – Thermomètre auriculaire EasyTherm</i>	C - 1
IT	<i>Istruzioni per l'uso – Termometro auricolare EasyTherm</i>	D - 1
NL	<i>Bedieningshandleiding – EasyTherm oorthermometer</i>	E - 1
ES	<i>Manual de instrucciones – Termómetro de oído EasyTherm</i>	F - 1
PT	<i>Manual de instruções – Termómetro auricular EasyTherm</i>	G - 1
PL	<i>Instrukcja obsługi – Termometr douzny EasyTherm</i>	H - 1
TR	<i>Kullanım kılavuzu – EasyTherm kulak termometresi</i>	I - 1
RU	<i>Инструкция по эксплуатации – ушной термометр EasyTherm</i>	J - 1
DA	<i>Betjeningsvejledning – luftfugtigheds- & temperaturmåler</i>	K - 1
FI	<i>Käyttöohje – EasyTherm korvakuumemittari</i>	L - 1
NO	<i>Bruksanvisning – EasyTherm øretermometer</i>	M - 1
SV	<i>Bruksanvisning – EasyTherm örontermometer</i>	N - 1

CE 0197

TROTEC® HealthCare GmbH • Grebbener Straße 7 • D-52525 Heinsberg
Tel.: +49 2452 962 - 530 • Fax: +49 2452 962 - 92 530
www.trotec-healthcare.de • E-Mail: info@trotec-healthcare.de

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b.) ASTM E1965-98(2009)

Standard Specification for Infrared Thermometers for Inter-mittent Determination of Patient Temperature

(ASTM: American Society for Testing and Materials)

Table of Contents

1. Model description.	B - 1
2. Safety instructions	B - 2
3. Unit depiction	B - 3
4. Functions	B - 3
5. Temperature measuring	B - 4
6. Display	B - 5
7. Instructions for use	B - 6
8. Changing temperature scales	B - 7
9. Troubleshooting	B - 8
10. Cleaning and maintenance	B - 9
11. Technical data	B-10


This infrared thermometer was built according to the latest technology. The device complies with the following standards: EN60825-1, EN61000 6-3 2001, EN 61000 6-1 2001 and fulfills the requirements of existing European and national guidelines. Please contact our service technicians on +49 2452 962 530 if you have any questions regarding electromagnetic tolerances or the measured values. Conformity has been proven. The pertinent declarations and documents are in the possession of the manufacturer.

1. Model Description

The digital ear thermometer you have purchased is a high-quality product that has been designed and developed in keeping with the highest and most modern technological standards and tested in compliance with international safety standards. The state-of-the-art technology used in the ear thermometer enables you to carry out precise and reliable measurements without any impairment to hearing during measuring. The thermometer automatically carries out an independent calibration check after powering on to ensure that you have an exact and reliable temperature measurement every single time.

Your digital ear thermometer is for private use only. It is designed for carrying out temperature measurements in intervals in order to provide you with a means of performing efficient and effective health checks. It is suitable for measuring the temperature of people of all ages.

Please read these Instructions carefully before using your ear thermometer and be sure to keep them in a safe and accessible place.

Electronic equipment may not be treated as domestic waste. In accordance with the EU (European Union) Directives on Waste Electrical and Electronic Equipment (WEEE), electronic equipment must be disposed of in a professional manner in accordance with Directive 2002/96EU of the EUROPEAN COUNCIL AND PARLIAMENT of 27 January 2003 with regard to old electrical and

electronic equipment. Please dispose of this appliance in a manner appropriate to the relevant legal requirements at the end of its product life.

2. Safety Instructions

- You should see a doctor if you have a temperature especially a high temperature over a longer period ! Compare the temperature you have taken using the measuring method you selected with the values below.

The temperatures you have taken with the ear thermometer are not a substitute for seeing a doctor.


- Your digital ear thermometer is not water resistant! Do not dip the thermometer in any liquid without first pulling over the protective sheath.
- Your ear thermometer may only be used for the purpose for which it was intended. Safety instructions must be adhered to at all times, especially where young children or babies are involved.
- Protect your digital thermometer against the sun and store in a dry and dust-free environment between 10° and 40° C (50° and 104°F).
- Do not use the thermometer if the measuring tip or the unit itself display even the slightest sign of damage. Do not on any account attempt to repair the thermometer yourself!

Contact a specialist immediately or call the following service number: +49 2452 962 - 530.

- Your thermometer is made of high-quality components which allow you to carry out accurate and reliable measurements! Handle with care! Do not drop!
- A collection of cerumen, commonly known as earwax, in the ear canal may lead to false readings. Please make sure that the ear canal of the person whose temperature you wish to measure is clean and free of cerumen.
- Protect against knocks and bumps. Do not try to twist or bend the thermometer or the measuring tip.
- In compliance with standards regarding safety-related checks and metrological inspections, it is recommended that the measuring device be subjected to a regular visual inspection of the housing components for any signs of damage and that the function of the keys and the measuring accuracy (see operating instructions) be checked once a year to ensure their correctness. The metrological inspection may either be carried out by the manufacturer or a specially-named or qualified testing centre.

3. Unit Depiction

1. Sensor
2. Start key
3. Display
4. On/Off key
5. Battery compartment cover

- = Identification number of the notified body
- = Observe accompanying documents!
- = for body temperature measurements
- = Classification for the disposal of electronic goods
- = Battery voltage 3V

4. Functions

Multiple Application Possibilities

You digital ear thermometer has an exceptionally broad measuring range (0° to 100° C = 32° to 212° F). This not only makes it especially suitable for measuring body temperature, but also for carrying out temperature measurements in the following cases:

- the surface temperature of milk in babies' bottles
- the surface temperature of the water in a baby bath
- the ambient air temperature

⚠ The surface temperature of liquids differs to the temperature at the core; the air temperature is only measured at exactly the spots where the measurements are carried out and is therefore not necessarily representative for the ambient air temperature.

The innovative infrared technology allows you to carry out measurements in seconds!

Accurate and reliable

The unique construction, the high-quality infrared sensor and the independent self-calibration function all ensure that the measurements are extremely accurate and reliable every single time.

Gentle and easy to use

The ear thermometer's ergonomic shape makes it especially easy to use. And because it is so easy to use one has no trouble measuring the temperature of sleeping children without waking them. The fact that the children can carry on sleeping is especially beneficial because it adds to the recovery process.

Temperature measurements can be carried out quickly and efficiently. This is especially practical when measuring the temperature of small children or babies for whom rectal or oral measurements can have something menacing about them.

Auto - Memory Funktion

The display shows the last measured temperature value for a full two seconds after the thermometer has powered up. The status message „MEM“ appears on the display during this time.

Clean and hygienic

- There is no risk of anyone hurting themselves on glass or coming into contact with toxic mercury after breakages.
- Highest possible levels of safety when measuring the temperatures of small children and babies.
- The tip of the thermometer can be cleaned quickly and easily using a cotton wool cloth that has previously been moistened with alcohol. The thermometer is now hygienically clean and can be used to measure the temperature of other family members.



Fever alarm

10 short alarm tones indicate that the child or adult whose temperature is being measured probably has fever.

5. Temperature Measurements

The digital infrared thermometer measures the infrared energy which is emitted from the ear drum and the surrounding tissue. This energy is collected through the lens and converted into a temperature value. Temperature measurements which are taken at the ear drum are the most precise of all. Temperature measurements which are taken exclusively from the surrounding tissue can result in lower temperature values. This in turn can lead to a fever being wrongly diagnosed.

Using The Thermometer Correctly In Order To Prevent Wrong or Inaccurate Readings

- Place a clean, new protective cap over the measuring tip before carrying out any measurements in the ear canal.

1. Take out a protective cap and place it in the holding fixture of the applicator included in the scope of delivery.

2. Press the measuring tip of the thermometer through the holding fixture of the applicator so that the protective cap fully encloses the measuring tip of the ear thermometer.

3. Remove the cap after use and replace with a new cap before carrying out further temperature especially when carrying out measurements on different persons.

- Press the on/off key to power up the thermometer. The thermometer emits a tone and the scale with the temperature units starts to blink. Pull the ear in question gently towards the back of the head and then upwards to straighten the ear canal. Place the sensor firmly (but not too firmly!) in the ear canal and press the start key. The sensor must remain in the ear canal until a further tone is emitted. This tone indicates that the temperature measurement is completed.
- Your digital ear thermometer has been clinically tested. The measurement results it provides are proven as long as it is used according to the operating instructions.

6. Display

Display	Meaning	Description
	All the symbols are shown.	Press the on/off key to power up. All the symbols appear for 2 seconds.
	Memory function	The last measurement is shown on the display for 2 seconds.
	Ready to use	The thermometer is ready to use. The temperature measurement can be carried out. The °C or °F symbol starts to blink.
	Measurement completed	The measured value is shown on the display. The backlit display lights up and then goes off and the °C or the °F symbol starts to blink. Temperature measurements can be resumed.
	The temperature was not measured inside the ear	The measured value is not within the measuring range of 32.0° to 42.2°C (89.6° to 108.0°F)
	Low battery	A blinking battery symbol indicates that the battery needs to be replaced.

7. Instructions For Use

Display	Meaning
	Press the on/off key. All the symbols appear for a period of 2 seconds on the display.
	The last measurement appears on the display together with the status message „MEM“ for a period of 2 seconds.
	The device is ready to use when the °C or °F symbol blinks on the display and the thermometer emits a tone.
	Pull the ear gently towards the back of the head and then upwards to straighten the ear canal in order to guarantee that the sensor has free passage to the ear drum. When measuring the temperature of infants under one year (babies), pull the ear gently towards the back of the head. In case of children, adolescents and adults, pull the ear gently backwards and then upwards.
	Carefully insert the sensor into the ear while the ear is gently being pulled towards the back of the head. Make sure that the sensor sits firmly (but not too firmly!) in the ear and press the start key. Release the start key and wait for the tone. This tone indicates that the temperature measurement has been completed.
	Remove the thermometer from the ear. The measured value is shown on the display.

- The thermometer emits 10 alarm tones when the measured temperature is higher than 37.5°C (99.5°F): This indicates that the person whose temperature has been measured may be suffering from fever. Wait at least 30 seconds between any two measurements if you wish to carry out 3 to 5 measurements so as to increase the accuracy of the measurements.
- A collection of cerumen („ear wax“) on the measuring tip of the thermometer can quickly lead to a wrong measurement or a cross-infection among those being examined. It is therefore essential that the measuring tip is cleaned before the next measurement is carried out. For more information, please see „Cleaning and Maintenance“.
- Wait at least 5 minutes before taking the next temperature when the measuring tip has been cleaned with alcohol. This gives the thermometer enough time to once again attain the reference value.
- The following method has proved to be especially successful when measuring the temperature of small children. Place the child on its back and turn the child’s head gently to one side so that the ear is facing upwards.
- When measuring the temperature of children or adults, we recommend that the temperature is taken from an angle behind the person’s ear.
- Always measure the temperature at the same ear. There are sometimes temperature differences between the right ear and the left ear.
- Please wait several minutes before taking the temperature of a person who has just fallen asleep.
- We recommend that you carry out three measurements in the following situations before selecting the highest value:
 - Newly-born babies during the first 100 days
 - Infants under 3 years with a compromised immune system when the onset of fever is critical
 - When the thermometer is being used for the first time and the person using it is unsure how it works and therefore needs to familiarize him or herself with it first.

8. Changing the Units of Temperature

Your digital ear thermometer can display temperatures both in degrees Celsius and Fahrenheit. To change from Celsius to Fahrenheit , simply power off the thermometer and press and hold the start key for 5 seconds. The selected temperature unit blinks on the display when the key is released. Press the start key again while the units are blinking to change from Celsius to Fahrenheit or vice versa. Wait another 5 seconds after the new units are displayed. The thermometer then automatically adjusts the temperature scale to the new units.

9. Troubleshooting

Display	Meaning	Possible Cause / Solution
	The temperature measured is too high	The symbol (H) appears in connection with (🔋) on the display when the temperature is lower than 100°C or 212°F.
	The temperature measured is too low	The symbol (L) appears in connection with (🔋) on the display when the temperature is lower than 0°C or 32°F.
	The operating temperature is too high	The symbol (H) appears on the display when the operating temperature is higher than 40°C or 104°F.
	The operating temperature is too low	The symbol (L) appears on the display when the operational temperature is lower than 10°C or 50°F.
	Error message	There is a system error! Check the battery and replace if necessary. Contact our Service Centre if the error cannot be rectified.
	Battery empty / Empty display	The batteries must be replaced immediately when the battery symbol is the only symbol on the display or the display is empty.

10. Cleaning and Maintenance

Cleaning and disinfection

Please check the housing of your infrared thermometer at regular intervals for any visible signs of damage or soiling and please make sure that the keys are in good working order. Please ensure that the sensor tip is clean and free of dirt or any other substances. Your infrared thermometer should be checked once a year using a suitable means of testing to check the measuring accuracy (see operating instructions).

Please use a 70% isopropyl alcohol solution and a moist swab or soft cloth to remove any dirt or other signs of fouling BEFORE taking the temperature. Please make sure that there are no remains of the cleaning solution on the thermometer when cleaning/disinfection has been completed.

⚠ Please ensure that no liquids are allowed to enter the thermometer casing. Do not use caustic cleaning agents, scouring powder or benzene to clean. Do not immerse in water or any other household detergent or cleaning agent.

Storage

Please store the thermometer at room temperature in its original casing or a protective sheath when not used over a longer period.

Please ensure that the display and lens of the measuring tip are not subjected to scratches.

Remove the batteries from the battery compartment when the device is not used or needed over a longer period. This prevents the batteries from causing any damage to the device in case of a leak.

Changing the battery

- The batteries are empty and need to be replaced when the battery symbol on the display starts to blink.
- Open the battery compartment and change the batteries. Remove the battery compartment cover by applying a little pressure to the rippled arrow symbol on the back of the thermometer and slide the cover off.
- Insert the batteries carefully and take care to make sure that the poles are properly aligned (+)/(-). Batteries which have not been inserted properly can damage your device. This can have an effect on the guarantee.
- Do not use rechargeable batteries. Use high-quality normal batteries only.

Lifetime

The infrared thermometer has been developed to meet the highest demands in daily and professional use. The lifetime is around 6 Years.

11. Technical Data

Mesasuring range..... Body temperature: 32,0° to 42,9°C
(89,6° to 109,2°F)

Extended temperature range: 0,0° to 100,0°C
(32,0° to 212,0°F)

Accuracy (lab conditions) ±0,2°C at 32,0° to 42,9°C
±0,4°F at 89,6° to 109,2°F
±0,1°C at 0,0° to 31,9°C; 43,0° to 100,0°C
±0,2°F at 32,0° to 89,5°F; 109,3° to 212,0°F;

Display..... LCD with increments of 0.1°C (0.1°F)

Acoustics..... The thermometer emits a short
tone when it is ready to measure.
A long tone is emitted when the measurement
is completed. Three short tones
indicate an error or a system error.

Memory function..... The last measured value
is automatically shown on the display.

Backlit display (optional)..... The backlit display
lights up for 4 seconds after powering on.
The backlit display goes off
approx. 5 seconds after the
measurement has been carried out.

Operating temperature.....10°C to 40°C (50°F to 104°F)

Storage temperature.....-25°C to 55°C (50°F to 104°F)

Auto power off..... after approx. 1 minute

Battery..... 2 x AAA batteries

Dimensions L x W x H.....155 x 33 x 45 mm
(6.1 x 1.3 x 1.77

ins)

Weight..... 87g (with battery); 68,5g (without battery)
(0.19/0.15

lbs)

Standards

a.) DIN EN 12470-5

Mecidinal thermometers – Section 5: Demands on
infrared ear thermometers (with maximum mechanism)

b.) ASTM E1965-98(2009)

Standard Specification for Infrared Thermometers for Inter-
mittent Determination of Patient Temperature

(ASTM: American Society for Testing and Materials)

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